


Spenco® PolySorb®

CROSSTRAINER

Replacement Insoles



Cat. No. 38-034

PERFORMANCE

Best combination of shock absorption and energy-return

Target Consumer:

Athletes looking for cushioning who participate in a variety of physical activities. Effective for use in all team sport competition and training. Consumers seeking additional cushioning in both the forefoot and heel.

Recommended For:

Increased shock absorption and heel strike protection. Reduced friction to help prevent blisters. Higher energy return for improved athletic performance.

Best For Use In:

Athletic, tennis, basketball and training shoes with removable insoles.

Components*:

PU, SpenCore™ Material, Fabric, Thickness at Heel: 13.21mm (1/2in), Thickness at Forefoot: 5.50mm (1/4in)

Instructions:

Remove the original insole from the shoe (some brands of shoes glue the insole to the shoe, but careful removal will not damage the shoe). Use the original insole as a guide to trace and trim the Spenco® PolySorb® Cross Trainer Replacement Insole.

Care:

Hand wash in warm water. Air dry.

Sizing Chart:

Size	0	1	2	3	4	5	6
U.S. Women's	3/4	5/6	7/8	9/10	11/12		
U.S. Men's			6/7	8/9	10/11	12/13	14/15
U.K.	2/3	4/5	5/6	7/8	9/10	11/12	13/14
European	32/34	36/38	38/40	40/42	42/44	44/46	46/48

Features & Benefits:

- 4-Way Stretch Fabric with Silpure Antimicrobial** helps prevent blisters while controlling odor.
- EVA Forefoot Cushion** provides superior cushioning and energy return for high impact sports.
- SpenCore™ Material** a 5/32" layer of this material absorbs shock and impact to provide intense comfort and cushioning.
- Lightweight POLYSORB™ Polyurethane Foam** provides arch and heel support.
- SpenCore™ Heel Plug** provides extra cushioning and shock absorption in the heel.

One-Year Unconditional Guarantee

Call toll-free: 1-800-877-3626
In Canada: 1-800-387-9538

www.spenco.com



*Measurements are approximate and vary by size.